

October 2007

Phone: 305-471-7670  
 Fax: 305-594-9035  
 sales@dloren-rasco.com



# SPLINT KITS

**Price**

S501-Half Arm	\$13.20
S502- Full Arm	\$15.10
S503-Half Leg	\$15.50
S504-Full Leg	\$18.80
S511-Foot & Ankle	\$14.40
S512-Hand & Wrist	\$11.90
S500-4 Splint Kit	\$59.20
S520 -6 Splint Kit	\$81.90
S530-4 Splint children	\$55.40

*Emergency First Aid pressure splints to immobilize injured limbs*

**Schuco Inflatable Air Splints**

Inflatable air splints provide immediate immobilization of injured limbs. Orally inflated via a screw air valve, the splint doubles as a self-contained tourniquet to stem the flow of blood and promote infection control. Transparent plastic and nylon zipper permit X-rays and constant observation of the injury. May be washed and sterilized by gas or cold solution for reuse.

**CLEAR PLASTIC ■ ORALLY INFLATED**

**DIRECTIONS FOR USING FIRST AID SPLINTS**

**INSERTION:**

- Select proper size splint and remove from protective bag (save bag).
- Open zip and carefully place splint around injured limb. Do not place zip over wound.
- To avoid splint punctures place gauze over compound fractures.

**INFLATION:**

- Close zip.
- Open valve. Inflate by mouth only, using full lung pressure. Close valve.

**CHECK PRESSURE:** If splint becomes too hard due to limb swelling, reduce pressure by releasing air. (see caution below.)

**REMOVAL:** Remove by opening valve and deflating splint. Open zip and remove splint.

**MAINTENANCE:** May be washed and cold sterilized for re-use. Light powdering with baby talc will extend life of splint.

**CAUTIONS:**

- Do not place zip over wound
- Inflate by mouth
- Inflate splint to point where your thumb will cause only a slight dent if pressed against the splint. Note that an injured limb may swell causing the splint to tighten and restrict blood circulation. Reduce air pressure so that circulation is restored while limb remains immobilized.

	SC 501 Half Arm
	SC 502 Full Arm
	SC 503 Half Leg
	SC 504 Full Leg
	SC 511 Foot & Ankle
	SC 512 Hand & Wrist
Cat No. SC 500B 4 Splint Kit (Adult) (Consists of: 501, 502, 503, 504, 505B & 506)	Splint Kit also available for Children
Cat No. SC 520B 6 Splint Kit (Adult) (Consists of: 501, 502, 503, 504, 511, 512, 505B, 506)	

Special points of interest:

**INDICATIONS**

Temporary Emergency Use for Fractures, Dislocations, Bruises, Sprains, Cuts, Burns and Related Injuries

**ADVANTAGES**

- Immobilizes injured limbs permitting safe transportation
- Stops bleeding, eases pain and prevents further injury and shock
- X-Ray transparent

D'Ioren-Rasco Corporation  
 10530 NW 26 Street, Suite F102  
 Doral, FL 33172 USA